

# Confident Quantification of Steroids: Analysis in Human Plasma or Serum by Liquid Chromatography-Tandem Mass Spectrometry for Clinical Research

Claudio De Nardi, Thermo Fisher Scientific, Dreieich, Germany; Sergio Indelicato, Thermo Fisher Scientific, Les Ulis, France; Maura Brambilla, Ospedale di Desio, Desio, Italy; Chiara Fania, Università degli Studi Milano-Bicocca, Milano, Italy

## ABSTRACT

**Purpose:** Development and implementation of a robust, reliable, sensitive analytical method for quantification of eight steroids in human plasma or serum using the Thermo Scientific™ TSQ Quantis™ triple quadrupole mass spectrometer.

**Methods:** Method performance was evaluated using the ClinMass® LC-MS/MS Complete Kit from RECIPE Chemicals + Instruments GmbH (Munich, Germany) in terms of linearity of response within the calibration ranges, carry over, accuracy and intra-assay precision for each analyte. Samples were extracted by protein precipitation followed by LC separation on a Thermo Scientific™ Transcend™ II LX-2 system using mobile phases and analytical column provided by RECIPE. Total runtime was 7.0 minutes. Analytes and internal standards were detected by selected reaction monitoring (SRM) on a TSQ Quantis triple quadrupole mass spectrometer with heated electrospray ionization operated in positive mode. Data were acquired and processed using Thermo Scientific™ TraceFinder™ 4.1 software.

**Results:** The method proved to be linear in the calibration ranges covered by the calibrators. It was not possible to detect the lowest calibrator for 17-hydroxyprogesterone, but the method was linear in the range calibrator 2 – calibrator 6. The maximum registered carry over was 0.08% with the exception of androstenedione, which showed a value of 0.36%. The data demonstrated outstanding accuracy of the method with the percentage bias between nominal and average back-calculated concentration for the used control samples ranging between -4.4% and 5.3%. The %CV for intra-assay precision was always below 5.7% for all the analytes.

## INTRODUCTION

The endocrine system controls signaling pathways that direct critical physiological processes such as development, reproduction, and metabolism. A major group of signaling molecules of the endocrine system is the steroid hormone family. Steroid hormones are defined as endogenous hormones synthesized from a cholesterol backbone and subtle changes in steroid structure can lead to drastic differences in physiological function. The four classes of steroid hormones (estrogens, androgens, progestogens, and corticosteroids) are often grouped together, however, most steroid hormones are regulated at different concentrations by different species. Quantifying and analyzing only one class of steroids often lead to misleading information and hence, there is a strong need for a comprehensive steroid hormone profile that include analysis and quantitation. While several technologies exist for quantitation of analytes in biological matrices, owing to its portability, ease-of-use, selectivity and specificity, liquid chromatography (LC) coupled to triple quadrupole mass spectrometers (MS/MS) has gained widespread popularity.

In this report, a robust, reliable LC-MS/MS method for clinical research is developed for the quantification of eight steroids in human plasma or serum. The analytical method reported in here includes 11- and 21-deoxycortisol, 17- and 21-hydroxyprogesterone, androstenedione, cortisol, dehydroepiandrosterone sulphate (DHEAS) and testosterone.

## MATERIALS AND METHODS

### Target analytes

The analytes and corresponding internal standards and concentration ranges covered by the calibrators used are reported in Table 1.

Table 1. Concentration ranges covered by the calibrators

Analyte	Internal Standard	Concentration Range (ng/mL)
Androstenedione	d7-Androstenedione	0.117 – 12.7
Cortisol	d4-Cortisol	2.75 – 325
DHEAS	d5-DHEAS	43.9 – 4851
11-Deoxycortisol	d5-11-Deoxycortisol	0.100 – 12.5
21-Deoxycortisol	d5-11-Deoxycortisol	0.117 – 14.7
17-Hydroxyprogesterone	d8-17-Hydroxyprogesterone	0.109 – 12.5
21-Hydroxyprogesterone	d8-17-Hydroxyprogesterone	0.196 – 24.4
Testosterone	d3-Testosterone	0.109 – 11.0

### Sample preparation

Reagents included seven calibrators (including blank) and three controls from RECIPE, as well as six isotopically-labelled internal standards for the quantification. Samples of 50 µL of plasma or serum were protein precipitated using 75 µL of precipitating solution containing the internal standards. Precipitated samples were vortex-mixed and centrifuged, and the supernatant was transferred to a clean plate or vial.

### Liquid chromatography

The LC separation was achieved using mobile phases and analytical column provided by RECIPE. Details of the analytical method are reported in Table 2. Total runtime was 7.0 minutes.

### Mass spectrometry

Analytes and internal standards were detected by selected reaction monitoring (SRM) on a TSQ Quantis triple quadrupole mass spectrometer with heated electrospray ionization operated in positive mode. Extracted samples were injected onto a Transcend II LX-2 system connected to TSQ Quantis followed by detection using six isotopically-labeled internal standards. A summary of the MS conditions is reported in Table 3. Two SRM transitions for each analyte were included in the acquisition method for quantification and confirmation, respectively.

Table 2. LC method description

	Time (min)	Flow Rate (mL/min)	A (%)	B (%)
Gradient Profile	0.00	0.5	73	27
	4.50	0.5	43	57
	4.51	0.5	0	100
	4.64	0.5	0	100
	4.75	0.8	0	100
	5.35	0.8	0	100
	5.36	0.8	73	27
	5.65	0.8	73	27
	5.66	0.5	73	27
	7.00	0.5	73	27
Injection Volume	40 µL	Column Temperature	40 °C	

Table 3. MS settings

Source type	Heated electrospray ionization
Vaporizer temperature	350 °C
Capillary temperature	325 °C
Spray voltage	3500 V (positive mode)
Sheath gas	50 AU
Sweep gas	1 AU
Auxiliary gas	10 AU
Data acquisition mode	Selected-reaction monitoring
Collision gas pressure	1.5 mTorr
Cycle time	0.300 s
Q1 mass resolution	0.7 (FWMH)
Q3 mass resolution	1.2 (FWMH)

### Method evaluation

Method performance was evaluated using the ClinMass® LC-MS/MS Complete Kit from RECIPE Chemicals + Instruments GmbH (Munich, Germany) in terms of linearity of response within the calibration ranges, carry over, accuracy and intra-assay precision for each analyte. Carry over was calculated in terms of percentage ratio between peak area of the highest calibrator and a blank sample injected just after it. Analytical accuracy was evaluated in terms of percentage bias between nominal and average back-calculated concentrations using quality control samples at three different levels provided by RECIPE (MS12083 batch #1226) prepared and analyzed in replicates of four. Intra-assay precision was evaluated on the same run (control samples at three levels, replicates of four) in terms of percentage coefficient of variation (%CV).

### Data analysis

Data were acquired and processed using TraceFinder 4.1 software.

## RESULTS

The method proved to be linear in the calibration ranges covered by the calibrators. It was not possible to detect the lowest calibrator for 17-hydroxyprogesterone but the method was linear in the range calibrator 2 – calibrator 6. Representative chromatograms for the lowest calibrator for testosterone, 21-deoxycortisol, and their internal standards are reported in Figure 1. Representative calibration curves for the same analytes are reported in Figure 2.

Figure 1. Representative chromatograms of the lowest calibrator for (a) testosterone, (b) d3-testosterone, (c) 21-deoxycortisol and (d) d5-11-deoxycortisol





Figure 2. Representative calibration curves for (a) testosterone and (b) 21-deoxycortisol



The maximum registered carry over was 0.08% with the exception of androstenedione which showed a value of 0.36%.

The data demonstrated outstanding accuracy of the method with the percentage bias between nominal and average back-calculated concentration for the used control samples ranging between -4.4% and 5.3%. The %CV for intra-assay precision was always below 5.7% for all the analytes. Results for accuracy and intra-assay precision reported in Table 4.

Table 4. Analytical accuracy and intra-assay precision results for control MS12083 batch #1226

Analyte	Level_I				Level_II				Level_III			
	Nominal Conc (ng/mL)	Average Calculated Conc (ng/mL)	Bias (%)	CV (%)	Nominal Conc (ng/mL)	Average Calculated Conc (ng/mL)	Bias (%)	CV (%)	Nominal Conc (ng/mL)	Average Calculated Conc (ng/mL)	Bias (%)	CV (%)
Androstenedione	0.219	0.223	1.9	3.7	0.604	0.604	0.0	2.5	3.21	3.29	2.6	0.8
Cortisol	5.53	5.61	1.5	1.0	15.4	16.0	3.7	0.9	84.0	88.5	5.1	1.1
DHEAS	86.4	89.0	2.9	5.2	235	241	2.3	4.0	1292	1345	4.0	2.4
11-Deoxycortisol	0.225	0.232	3.0	2.1	0.633	0.662	4.3	3.9	3.54	3.70	4.3	1.1
21-Deoxycortisol	0.248	0.256	3.0	2.1	0.692	0.699	1.0	3.9	3.76	3.97	5.3	0.5
17-Hydroxyprogesterone	N/A	N/A	N/A	N/A	0.599	0.574	-4.4	1.8	3.29	3.47	5.2	0.0
21-Hydroxyprogesterone	0.422	0.412	-2.5	4.4	1.18	1.18	-0.4	5.0	6.45	6.80	5.1	1.5
Testosterone	0.198	0.196	-1.0	5.7	0.540	0.546	1.1	1.6	2.88	3.04	5.2	1.3

## CONCLUSIONS

A robust, reliable LC-MS/MS method for clinical research for the quantification of eight steroids in human plasma or serum was developed and implemented. The ClinMass LC-MS/MS Complete Kit from RECIPE ensured increased confidence in the result that was obtained. The method was analytically validated on a Transcend II system connected to a TSQ Quantis triple quadrupole mass spectrometer. The sample preparation procedure utilized in this method uses quick and simple offline protein precipitation with concomitant internal standard addition. The described method meets research laboratory requirements in terms of sensitivity, linearity of response, accuracy, and precision.

## ACKNOWLEDGEMENTS

We would like to acknowledge Katharina Kern from RECIPE for the support provided during the implementation of this application.

## TRADEMARKS/LICENSES

© 2018 Thermo Fisher Scientific Inc.

ClinMass is a registered trademark of RECIPE Chemicals + Instruments GmbH.

All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

## FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.