

Preliminary evaluation of the Thermo Scientific Cascadion SM LC-MS/MS random access clinical analyser in a fully automated routine clinical laboratory

Godwin K Tetteh, Sarah-Jayne Needham, Steven Alderson and Sally C Benton

Department of Blood Sciences, Berkshire & Surrey Pathology Services, Frimley Health NHS Foundation Trust

Introduction

- Liquid chromatography mass spectrometry (LCMS) is considered the gold standard method for analysis of a number of key analytes because of its enhanced analytical specificity relative to immunoassay.
- Despite degrees of automation within systems, LCMS is still a labour intensive technique for clinical laboratories requiring large amounts of staff time to carry out sample preparation before analysis.
- The Thermo Scientific™ Cascadion™ SM Clinical Analyser (Cascadion) Fig 1, is a fully automated LCMS random access analyser that has been recently launched. The system incorporates on-line sample preparation and LCMS technology in a closed analytical system.
- Preliminary assays planned for availability on the Cascadion are total vitamin D (reported as a total of D2 + D3), an immunosuppressant panel (tacrolimus, sirolimus, cyclosporine, everolimus) and testosterone. Currently the total vitamin D is CE-IVD compliant.

Aim

- To carry out an independent evaluation of the Cascadion and to assess its acceptability within a busy routine clinical laboratory with no prior LCMS experience
- To carry out an analytical validation of the vitamin D3 assay component only due to the rarity of D2 bearing samples.



Figure 1: Cascadion at Frimley Park Hospital

Method

- Two state registered biomedical scientists with no previous experience of LCMS were selected as key operators. They received 3 days off-site training by Thermo Fisher Scientific and 2 days of on-site refresher training.
- Upon installation of the Cascadion in the laboratory the two scientists were tasked with running routine serum samples through the analyser for vitamin D analysis.
- User acceptability was defined as the ability to independently carry out this analysis as well as carry out routine maintenance.
- Vitamin D3 validation
 - Vitamin D stripped serum was run 20 times consecutively to assess limit of blank.
 - Carryover was assessed by analysing a high sample (D3 >301 nmol/L) in duplicate followed by a PBS sample in duplicate. This was repeated 10 times.
 - A high sample was diluted 1:2, 1:4, 1:8 and 1:16 to assess linearity.
 - Pooled serum samples were analysed 10 times consecutively to assess within batch precision and on 5 consecutive days to assess between batch imprecision.

Results

- The Cascadion is fully automated from sample preparation to result generation making it easy to operate.
- The two key operators had an automation background and no prior experience of LCMS. Following the 3 day off-site training and 2 day on-site refresher training staff were confidently able to run the Cascadion, carry out routine maintenance and perform basic trouble-shooting.
- The analyser has an automatic start-up. Routine daily maintenance takes a maximum of 10 minutes and routine monthly maintenance takes a maximum of 60 minutes
- Vitamin D3 validation
 - The vitamin D3 assay has an analytical range of 3.4 – 132 ng/mL (8.49 - 329.47 nmol/L).
 - No response was observed on a blank sample
 - No carryover was observed
 - Linearity of the assay was acceptable ($r^2 = 0.9995$)
 - Acceptable within batch and between batch precision was observed (tables 2a and 2b)

Sample	Mean ng/mL	SD	CV(%)
Low	19.13	0.57	2.96
Medium	87.15	3.94	4.53
High	131.16	4.01	3.06

Table 2a: Between batch precision

Sample	Mean ng/mL	SD	CV(%)
Low	18.60	0.95	5.10
Medium	87.93	3.22	3.66
High	128.10	1.93	1.51

Table 2b: Within batch precision

Conclusion

- In this study we have demonstrated that state registered scientists with no experience of LCMS are able to confidently run the Cascadion analyser with minimal training.
- Preliminary data for the D3 assay has demonstrated acceptable performance.
- A sample comparison for vitamin D3 and validation of the vitamin D2 assay is on-going.
- Further work is continuing to fully evaluate the Cascadion system including an assessment of sample throughput and walkaway time for 24-7 working.

Product is IVD/CE marked but not 510(k) –cleared and not yet available for sale in the U.S. Availability of product in each country depends on local regulatory marketing authorization status.

Acknowledgements

We are grateful to the entire Thermo Fisher team and the Blood Sciences team at Frimley Park for all their support with this work.